



always onestepahead



Children's 2018

Made to Order

Paediatric Therapeutic Footwear

Piedro[®]

is one of the world's leading
manufacturers of high-quality
Therapeutic Footwear.

Welcome to the Piedro® Made to Order Range of Paediatric Therapeutic Footwear.

Piedro have been manufacturing footwear for over 60 years. In that time the name Piedro has become the benchmark for quality, fit, durability and style in paediatric footwear. Our latest range reflects superb cosmesis and high street styling. We have used colourful leathers, fun prints and modern designs to make a contemporary range of great children's footwear. We have a unique collection of mid tops and high tops plus our trainer range which offer a modern option for casual wear. Plus, we have a selection of footwear which is ideal for school wear.



Piedro Footwear is made to four fabrication types:

Rehabilitation

Stability

AFO

Reverse Last (Anti-Varus)

See pages 8 - 9 to view how each type is constructed.

Most Piedro footwear designs come in all 4 different fabrication types, with some exceptions, which are clearly indicated on the pages concerned.

All Piedro footwear is made to exacting standards using the highest levels of craftsmanship, quality materials and cutting-edge technology.

Our leathers are specially selected for durability and all our footwear has a superb calfskin lining.

We are renowned for the quality of our products, and therefore, we guarantee our footwear for 6 months (excluding normal wear and tear).

NEW. We are delighted to announce our Made to Order range of **Piedro®** footwear.

If you foresee a problem and want to change the fit of a standard Piedro boot, **you can change it.**

You can **order split sizes** or **increase the volume and shape of our Rehabilitation, Stability, AFO and Reverse Last** footwear.

You can also change:

- the sole unit you require
- fastenings
- colours
- the style of toe puff
- the strength of heel counter
- you can shape heel counters around a Haglund's exostosis or prominent malleolus
- ask for extra ankle padding
- or request a cut out to be included in the ankle stiffener

Plus many other variations on our standard models.

Please advise us of where on the last and the size of last additions you require, **see page 85**, here you can also view other alterations which may be made to our footwear.

If what you want isn't listed then please talk to our customer services team who will be happy to help.

Piedro Customer Services
0845 617 1976

Delivered to you within 15 - 21 working days.

(Footwear requiring fitting stage will be delivered finished within a further 15 working days).

All of our footwear can have orthotic adaptations applied to them, please see pages 81 - 83 for further instructions and codes.

All of our footwear can be fitted with our range of Functional Foot Orthoses which are made to our lasts and bottom patterns to ensure perfect fit, please see pages 81 - 82 for a full description of Functional Foot Orthoses available.



Please ensure a speedy process for despatching your order by stating clearly your requirements, footwear type, colour and footwear reference number, name of the shoe.

Piedro U.K. has a fully dedicated Customer Services Team based in Chesterfield who will always be pleased to assist with any of your enquiries. We also have a Business Development Manager who has many years of experience both clinically and technically who can answer any questions about fit, function and suitability of our range.

How to order the Made to Order footwear you require

The delivery time of Made to Order Footwear under normal circumstances is **15 - 21 working days**. Should you require a fitting stage **we will require a further 10 working days** to finish your order from receipt of the footwear into our offices.

DUE TO THE PERSONALIZED NATURE OF MADE TO ORDER FOOTWEAR, THEY ARE TREATED AS NON-RETURNABLE.

If you have any questions please ring us, we'll be happy to help.

0845 617 1976



Order process

- 1 Select the fabrication type of footwear you require

Rehabilitation

Stability

AFO

Reverse Last (Anti-Varus)

from our **Made to Order** range

- 2 Select your size and width
- 3 Select the style of footwear you require
- 4 Select the Name and Reference number of the style you require

Name

Graff

Reference number

2034.9957

- 5 Tell us where you want additions/changes made to the standard Piedro footwear.
(see chart on page 85 for guidance if required)
- 6 Choose your Orthotic adaptations as required
- 7 Contact us with your requirements



customerservices@piedro-uk.co.uk



0845 617 1976



0845 617 1977



piedro-uk.co.uk

always one step ahead



PIEDRO





Made to Order Contents

8 - 9	Construction Styles
10 - 11	How to Measure Piedro Footwear
12 - 16	Footwear Size and Width Conversion Charts
17 - 18	Footwear Adaptions and Repairs

Made to Order Collection

20 - 23	Firststeps Sizes 15 - 21
24 - 47	Stitchdown Sizes 18 - 42
48 - 52	Trainers Sizes 24 - 42
54 - 60	High Tops & Mid Tops Sizes 24 - 42
62 - 69	Sandals Sizes 24 - 46
70 - 76	Only AFO Collection Sizes 24 - 46

Insoles, Accessories and Footwear Care

80	Heel Cups
81 - 82	BioFit Functional Inlays
83	Respro Standard and Arch Filled Inlays
84	Accessories
85	Order Form
86	Order Information

Construction styles

Rehabilitation



Features:

- Made on extra depth lasts
- Extended medial and lateral heel stiffeners to provide extra foot support
- Removable, anatomically shaped footbeds
- Adaptable
- Extended heel stiffeners

Indications for use:

- Children who need footwear with extra support
- Neuro-physiological problems
- Varus or valgus instability
- Painful or pressure sensitive feet
- To accommodate insoles

1. Soft padding around the collar for a better fit
2. Extended stiffeners for better medial and lateral stability
4. Full leather insole
5. Easily removable and adaptable (hard-wearing) sole-unit
6. Sufficient toe-spring for easy roll-off
7. Removable profiled cork inlay, inside shoe
8. High toe-box and strong toe-cap for protection of the toes
9. Extended facings for maximum accessibility
10. Extra padding around the ankle for more comfort
11. Soft tongue padding to mould the shoe round the foot

Stability



Features:

- Fixates the subtalar joints and stabilises the talocrural joint
- Made on extra depth lasts
- High and extended heel stiffeners
- Removable, anatomical footbeds
- Adaptable

Indications for use:

- Children who need footwear with extra support
- Feet with extreme valgus or varus instability
- Severe Cerebral Palsy (CP)
- Partial or complete paralysis of the malleolus area
- Limited dorsi-flexion
- Deterioration of the talocrural joint
- Problems with ligaments in the lower limbs
- Painful or pressure sensitive feet
- To accommodate insoles

1. Soft padding around the collar for a better fit
2. A high and extra long stiffener to support the ankle
3. Full leather insole
4. Easily removable and adaptable (hard-wearing) sole-unit
5. Sufficient toe-spring for easy roll-off
6. Removable profiled cork inlay, inside shoe
7. High toe-box and strong toe-cap for protection of the toes
8. Extended facings for maximum accessibility
9. Extra padding around the ankle for more comfort
10. Soft tongue padding to mould the shoe round the foot

AFO

Reverse Last



Features:

- Made on specialist footwear lasts to ensure the correct position of the lower limb is at 90°
- Additional space provided at the front and heel
- Can accommodate stock or custom Ankle Foot Orthosis (AFO)
- Adaptable

Indications for use:

- To accommodate an Ankle Foot Orthosis (AFO)
- To accommodate an inner-shoe such as a Rainier Boot or Arizona splint for example

1. Soft padding around the collar to fix the AFO
2. Special shaped extended heel stiffener moulded to fit the AFO
3. Full leather flat insole to allow the AFO to sit at a 90° angle
4. Easily removable and adaptable lightweight sole-unit with anti-pronation function
5. Hard-wearing removable anti-slip full rubber sole layer
6. Sufficient toe-spring for easy roll-off
7. High toe-box and strong toe-cap for protection of the toes
8. Extended facings for maximum accessibility
9. Soft tongue padding to mould the shoe round the foot
10. AFO



Features:

- Specialist last shape
- Extended medial stiffener to prevent the foot progressing further into adduction
- Made on extra depth lasts
- Removable anatomical footbeds
- Adaptable

Indications for use:

- Clubfoot or to stabilise the post-operative clubfoot
- Extreme varus and / or adduction of the foot
- To correct the peak position in the upper tarsal joint (*pes equino-varus*)
- To correct varus tilting in the lower tarsal joint and further metatarsal area (*pes varus*)
- To stimulate abduction of the mid and forefoot
- To improve toe-off

1. Soft padding around the collar for a better fit
2. A special extended stiffener with curved wall to correct the foot shape
3. Full leather insole
4. Easily removable and adaptable (hard-wearing) sole-unit
5. Sufficient toe-spring for easy roll-off
6. Removable profiled cork inlay, inside shoe
7. High toe-box and strong toe-cap for protection of the toes
8. Extended facings for maximum accessibility
9. Extra padding around the ankle for more comfort
10. Soft tongue padding to mould the shoe round the foot

Measure your shoe size and foot width

Piedro® footwear is available in 32 sizes, from 15 to 46, with most shoes available in up to 8 width fittings. It is very easy to work out exactly what size and width you need using our conversion charts table on pages 12 to 16.

How to find your Shoe Size:

How to find your Shoe Size: Look at the four diagrams (Fig. A, Fig. B, Fig. C & Fig. D) to see where you need to measure for foot circumference, width and size. Take these four measurements (as shown) in millimetres and make a note of them for later. Take your size measurement using the charts on pages 12 - 14, scroll down the left hand 'Foot Size' column until you see the closest measurement to the one you have taken.

Next, work your way across the page until you find the closest width size to match your width measurement. Follow this all the way to the top column and this will give you the width required. Now, you have all of the measurements you need to make your Piedro® Footwear selection.

If you need any help or advice with Piedro® footwear sizing, please contact our Customer Services Department on **0845 617 1976**. They will be glad to help you.

Using the Piedro® Foot Measure (Code **99001**) makes the need to account for growth or weight bearing variances unnecessary.

Split Sizes can be supplied from stock, surcharge applies.



Code: **99001**
Piedro® Foot Measure
Measure widths 1½ to 8½




Fig A. **Joint circumference**




Fig D. **Short Heel circumference**




Fig B. **Joint Width measurement**




Fig C. **Size measurement**

How to measure for Piedro® Footwear



Using the Piedro® Calibrated Foot Measure

1. Position the foot on the measure non-weight bearing and move the slider up to the longest toe
2. Read the shoe size required, in green by the arrow
3. Move the tape slider to the widest part of the foot and measure the joint circumference around this point
4. Next, reposition the foot sideways, then, move the slider up to the widest part of the foot and measure the joint width in mm in the pink section
5. Compare these 3 measurements, size first, to the conversion chart and choose the closest shoe width to your measurements. For AFO Footwear, measure over the AFO and use the AFO conversion chart



Using the Piedro® Calibrated Foot Stick

1. Position the foot on the measure non-weight bearing and move the slider up to the longest toe
2. Read the shoe size required, in green by the arrow
3. Move the tape slider to the widest part of the foot and measure the joint circumference around this point
4. Next, reposition the foot sideways, then, move the slider up to the widest part of the foot and measure the joint width in mm in the pink section
5. Compare these 3 measurements, size first, to the conversion chart and choose the closest shoe width to your measurements. For AFO Footwear, measure over the AFO and use the AFO conversion chart



Using an alternative Non Calibrated Foot Stick

1. Position the foot on the stick non-weight bearing and move the slider up to the longest toe
2. Read the UK foot size and add 3 sizes (24mm). Move the slider down to this measurement and read the European shoe size required from the scale on the left. When measuring weight bearing, add only 2 sizes and follow the same process
3. Wrap a tape measure around to the widest part of the foot and measure the joint circumference around this point
4. Next, reposition the foot sideways, then, move the slider up to the widest part of the foot and measure the joint width in mm from the back of the stick
5. Compare these 3 measurements, size first, to the conversion chart and choose the closest shoe width to your measurements

Piedro® Footwear Size and Width Conversion Chart

Rehabilitation

Stability

Reverse Last

Shoe Size		Average Foot Size	Width 1.5						Width 2.5					
UK	EU		Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)
0.5	15	90	41.00	112.50	114.00	149.40	13.30	3.20	43.00	117.50	119.00	153.40	13.60	3.20
1	16	96	43.00	116.50	118.00	157.00	13.70	3.75	45.00	121.50	123.00	161.00	14.00	3.75
1.5	17	102	44.00	120.50	122.00	164.60	14.10	4.30	46.00	125.50	127.00	168.60	14.40	4.30
2.5	18	109	46.00	124.50	126.00	172.20	14.50	4.85	48.00	129.50	131.00	176.20	14.80	4.85
3	19	115	48.00	128.50	130.00	179.80	14.90	5.40	50.00	133.50	135.00	183.80	15.20	5.40
4	20	121	49.00	132.50	134.00	187.40	15.30	5.95	51.00	137.50	139.00	191.40	15.60	5.95
5	21	128	51.00	136.50	138.00	195.00	15.70	6.50	53.00	141.50	143.00	199.00	16.00	6.50
5.5	22	134	53.00	140.50	142.00	202.60	16.10	7.05	55.00	145.50	147.00	206.60	16.40	7.05
6	23	140	55.00	144.50	146.00	210.20	16.50	7.60	57.00	149.50	151.00	214.20	16.80	7.60
7	24	147	56.00	148.50	150.00	217.80	16.90	8.15	58.00	153.50	155.00	221.80	17.20	8.15
8	25	153	57.00	152.50	154.00	225.40	17.30	8.70	59.00	157.50	159.00	229.40	17.60	8.70
8.5	26	159	59.00	157.50	158.00	233.00	17.70	9.25	61.00	162.50	163.00	237.00	18.00	9.25
9.5	27	166	60.00	161.50	162.00	240.60	18.10	9.80	63.00	166.50	167.00	244.60	18.40	9.80
10	28	173	62.00	165.50	166.00	248.20	18.50	10.35	64.00	170.50	171.00	252.20	18.80	10.35
11	29	179	63.00	169.50	170.00	255.80	18.90	10.90	65.00	174.50	175.00	259.80	19.20	10.90
12	30	185	65.00	173.00	173.50	263.40	19.30	11.45	66.00	178.00	178.50	267.40	19.60	11.45
12.5	31	192	66.00	176.50	177.00	271.00	19.70	12.00	68.00	181.50	182.00	275.00	20.00	12.00
13.5	32	198	67.00	180.00	180.50	278.60	20.10	12.55	69.00	185.00	185.50	282.60	20.40	12.55
1	33	205	69.00	184.00	184.50	286.20	20.50	13.10	71.00	189.00	189.50	290.20	20.80	13.10
2	34	212	70.00	188.00	188.50	293.80	20.90	13.65	72.00	193.00	193.50	297.80	21.20	13.65
2.5	35	218	72.00	192.00	192.50	301.40	21.30	14.20	74.00	197.00	197.50	305.40	21.60	14.20
3.5	36	225	73.00	196.50	196.50	309.00	21.70	14.75	75.00	201.00	201.50	313.00	22.00	14.75
4	37	232	75.00	201.00	200.50	316.60	22.10	15.30	77.00	206.00	205.50	320.60	22.40	15.30
5	38	238	76.00	205.00	204.50	324.20	22.50	15.85	78.00	210.00	209.50	328.20	22.80	15.85
6	39	245	78.00	209.00	208.50	331.80	22.90	16.40	80.00	214.00	213.50	335.80	23.20	16.40
6.5	40	252	79.00	213.50	212.50	339.40	23.30	16.95	81.00	218.50	217.50	343.40	23.60	16.95
7.5	41	258	81.00	217.50	216.50	347.00	23.70	17.50	83.00	222.50	221.50	351.00	24.00	17.50
8	42	265	83.00	222.00	220.50	354.60	24.10	18.05	85.00	227.00	225.50	358.60	24.40	18.05
9	43	272	84.00	226.00	224.50	362.20	24.50	18.60	87.00	231.00	229.50	366.20	24.80	18.60
10	44	279	86.00	230.00	228.50	369.80	24.90	19.15	88.00	235.00	233.50	373.80	25.20	19.15
10.5	45	285	88.00	234.00	232.50	377.40	25.30	19.70	89.00	239.00	237.50	381.40	25.60	19.70
11.5	46	292	89.00	238.00	236.50	385.00	25.70	20.25	91.00	243.00	241.50	389.00	26.00	20.25

Rehabilitation

Stability

Reverse Last

Shoe Size		Average Foot Size	Width 3.5						Width 4.5					
UK	EU	Size (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)
0.5	15	90	45.00	122.50	124.00	157.40	13.90	3.20	47.00	127.50	129.00	161.40	14.20	3.20
1	16	96	47.00	126.50	128.00	165.00	14.30	3.75	48.00	131.50	133.00	169.00	14.60	3.75
1.5	17	102	48.00	130.50	132.00	172.60	14.70	4.30	50.00	135.50	137.00	176.60	15.00	4.30
2.5	18	109	50.00	134.50	136.00	180.20	15.10	4.85	52.00	139.50	141.00	184.20	15.40	4.85
3	19	115	52.00	138.50	140.00	187.80	15.50	5.40	54.00	143.50	145.00	191.80	15.80	5.40
4	20	121	53.00	142.50	144.00	195.40	15.90	5.95	55.00	147.50	149.00	199.40	16.20	5.95
5	21	128	55.00	146.50	148.00	203.00	16.30	6.50	57.00	151.50	153.00	207.00	16.60	6.50
5.5	22	134	57.00	150.50	152.00	210.60	16.70	7.05	58.00	155.50	157.00	214.60	17.00	7.05
6	23	140	58.00	154.50	156.00	218.20	17.10	7.60	60.00	159.50	161.00	222.20	17.40	7.60
7	24	147	59.00	158.50	160.00	225.80	17.50	8.15	61.00	163.50	165.00	229.80	17.80	8.15
8	25	153	61.00	162.50	164.00	233.40	17.90	8.70	63.00	167.50	169.00	237.40	18.20	8.70
8.5	26	159	63.00	167.50	168.00	241.00	18.30	9.25	65.00	172.50	173.00	245.00	18.60	9.25
9.5	27	166	65.00	171.50	172.00	248.60	18.70	9.80	66.00	176.50	177.00	252.60	19.00	9.80
10	28	173	66.00	175.50	176.00	256.20	19.10	10.35	68.00	180.50	181.00	260.20	19.40	10.35
11	29	179	67.00	179.50	180.00	263.80	19.50	10.90	69.00	184.50	185.00	267.80	19.80	10.90
12	30	185	69.00	183.00	183.50	271.40	19.90	11.45	71.00	188.00	188.50	275.40	20.20	11.45
12.5	31	192	70.00	186.50	187.00	279.00	20.30	12.00	72.00	191.50	192.00	283.00	20.60	12.00
13.5	32	198	71.00	190.00	190.50	286.60	20.70	12.55	73.00	195.00	195.50	290.60	21.00	12.55
1	33	205	73.00	194.00	194.50	294.20	21.10	13.10	75.00	199.00	199.50	298.20	21.40	13.10
2	34	212	74.00	198.00	198.50	301.80	21.50	13.65	76.00	203.00	203.50	305.80	21.80	13.65
2.5	35	218	76.00	202.00	202.50	309.40	21.90	14.20	78.00	207.00	207.50	313.40	22.20	14.20
3.5	36	225	77.00	206.00	206.50	317.00	22.30	14.75	79.00	211.00	211.50	321.00	22.60	14.75
4	37	232	79.00	211.00	210.50	324.60	22.70	15.30	81.00	216.00	215.50	328.60	23.00	15.30
5	38	238	80.00	215.00	214.50	332.20	23.10	15.85	82.00	220.00	219.50	336.20	23.40	15.85
6	39	245	82.00	219.00	218.50	339.80	23.50	16.40	84.00	224.00	223.50	343.80	23.80	16.40
6.5	40	252	84.00	223.50	222.50	347.40	23.90	16.95	86.00	228.50	227.50	351.40	24.20	16.95
7.5	41	258	85.00	227.50	226.50	355.00	24.30	17.50	87.00	232.50	231.50	359.00	24.60	17.50
8	42	265	87.00	232.00	230.50	362.60	24.70	18.05	89.00	237.00	235.50	366.60	25.00	18.05
9	43	272	89.00	236.00	234.50	370.20	25.10	18.60	90.00	241.00	239.50	374.20	25.40	18.60
10	44	279	90.00	240.00	238.50	377.80	25.50	19.15	91.00	245.00	243.50	381.80	25.80	19.15
10.5	45	285	91.00	244.00	242.50	385.40	25.90	19.70	93.00	249.00	247.50	389.40	26.20	19.70
11.5	46	292	93.00	248.00	246.50	393.00	26.30	20.25	95.00	253.00	251.50	397.00	26.60	20.25

Piedro® Footwear Size and Width Conversion Chart

Rehabilitation

Stability

Reverse Last

Shoe Size		Average Foot Size	Width 5.5						Width 6.5					
UK	EU		Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)
0.5	15	90	49.00	132.50	134.00	165.40	14.50	3.20	51.00	137.50	139.00	171.40	14.80	3.20
1	16	96	50.00	136.50	138.00	173.00	14.90	3.75	52.00	141.50	143.00	179.00	15.20	3.75
1.5	17	102	52.00	140.50	142.00	180.60	15.30	4.30	54.00	145.50	147.00	186.60	15.60	4.30
2.5	18	109	54.00	144.50	146.00	188.20	15.70	4.85	56.00	149.50	151.00	194.20	16.00	4.85
3	19	115	55.00	148.50	150.00	195.80	16.10	5.40	57.00	153.50	155.00	201.80	16.40	5.40
4	20	121	57.00	152.50	154.00	203.40	16.50	5.95	59.00	157.50	159.00	209.40	16.80	5.95
5	21	128	58.00	156.50	158.00	211.00	16.90	6.50	60.00	161.50	163.00	217.00	17.20	6.50
5.5	22	134	60.00	160.50	162.00	218.60	17.30	7.05	62.00	165.50	167.00	224.60	17.60	7.05
6	23	140	62.00	164.50	166.00	226.20	17.70	7.60	64.00	169.50	171.00	232.20	18.00	7.60
7	24	147	63.00	168.50	170.00	233.80	18.10	8.15	65.00	173.50	175.00	239.80	18.40	8.15
8	25	153	65.00	172.50	174.00	241.40	18.50	8.70	67.00	177.50	179.00	247.40	18.80	8.70
8.5	26	159	67.00	177.50	178.00	249.00	18.90	9.25	68.00	182.50	183.00	255.00	19.20	9.25
9.5	27	166	68.00	181.50	182.00	256.60	19.30	9.80	70.00	186.50	187.00	262.60	19.60	9.80
10	28	173	70.00	185.50	186.00	264.20	19.70	10.35	71.00	190.50	191.00	270.20	20.00	10.35
11	29	179	71.00	189.50	190.00	271.80	20.10	10.90	73.00	193.90	195.00	277.80	20.40	10.90
12	30	185	73.00	193.00	193.50	279.40	20.50	11.45	74.00	198.00	198.50	285.40	20.80	11.45
12.5	31	192	74.00	196.50	197.00	287.00	20.90	12.00	76.00	201.50	202.00	293.00	21.20	12.00
13.5	32	198	75.00	200.00	200.50	294.60	21.30	12.55	78.00	205.00	205.50	300.60	21.60	12.55
1	33	205	77.00	204.00	204.50	302.20	21.70	13.10	79.00	209.00	209.50	308.20	22.00	13.10
2	34	212	78.00	208.00	208.50	309.80	22.10	13.65	81.00	213.00	213.50	315.80	22.40	13.65
2.5	35	218	80.00	212.00	212.50	317.40	22.50	14.20	82.00	217.00	217.50	323.40	22.80	14.20
3.5	36	225	81.00	216.00	216.50	325.00	22.90	14.75	83.00	221.00	221.50	331.00	23.20	14.75
4	37	232	83.00	221.00	220.50	332.60	23.30	15.30	85.00	226.00	225.50	338.60	23.60	15.30
5	38	238	84.00	225.00	224.50	340.20	23.70	15.85	86.00	230.00	229.50	346.20	24.00	15.85
6	39	245	86.00	229.00	228.50	347.80	24.10	16.40	88.00	234.00	233.50	353.80	24.40	16.40
6.5	40	252	88.00	233.50	232.50	355.40	24.50	16.95	89.00	238.50	237.50	361.40	24.80	16.95
7.5	41	258	89.00	237.50	236.50	363.00	24.90	17.50	91.00	242.50	241.50	369.00	25.20	17.50
8	42	265	90.00	242.00	240.50	370.60	25.30	18.05	92.00	247.00	245.50	376.60	25.60	18.05
9	43	272	91.00	246.00	244.50	378.20	25.70	18.60	94.00	251.00	249.50	384.20	26.00	18.60
10	44	279	93.00	250.00	248.50	385.80	26.10	19.15	96.00	255.00	253.50	391.80	26.40	19.15
10.5	45	285	95.00	254.00	252.50	393.40	26.50	19.70	97.00	259.00	257.50	399.40	26.80	19.70
11.5	46	292	97.00	258.00	256.50	401.00	26.90	20.25	98.00	263.00	261.50	407.00	27.20	20.25

Rehabilitation

Stability

Reverse Last

Shoe Size		Average Foot Size	Width 7.5						Width 8.5					
UK	EU	Size (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)	Joint Width (mm)	Joint Circ. (mm)	Instep Circ. (mm)	Short Heel Circ. (mm)	Toe Box Height (mm)	Heel Height (mm)
0.5	15	90	53.00	142.50	144.00	175.40	15.10	3.20	55.00	147.50	149.00	179.40	15.40	3.20
1	16	96	54.00	146.50	148.00	183.00	15.50	3.75	56.00	151.50	153.00	187.00	15.80	3.75
1.5	17	102	56.00	150.50	152.00	190.60	15.90	4.30	58.00	155.50	157.00	194.60	16.20	4.30
2.5	18	109	57.00	154.50	156.00	198.20	16.30	4.85	59.00	159.50	161.00	202.20	16.60	4.85
3	19	115	59.00	158.50	160.00	205.80	16.70	5.40	61.00	163.50	165.00	209.80	17.00	5.40
4	20	121	60.00	162.50	164.00	213.40	17.10	5.95	62.00	167.50	169.00	217.40	17.40	5.95
5	21	128	62.00	166.50	168.00	221.00	17.50	6.50	64.00	171.50	173.00	225.00	17.80	6.50
5.5	22	134	64.00	170.50	172.00	228.60	17.90	7.05	66.00	175.50	177.00	232.60	18.20	7.05
6	23	140	65.00	174.50	176.00	236.20	18.30	7.60	67.00	179.50	181.00	240.20	18.60	7.60
7	24	147	67.00	178.50	180.00	243.80	18.70	8.15	69.00	183.50	185.00	247.80	19.00	8.15
8	25	153	68.00	182.50	184.00	251.40	19.10	8.70	70.00	187.50	189.00	255.40	19.40	8.70
8.5	26	159	70.00	187.50	188.00	259.00	19.50	9.25	72.00	192.50	193.00	263.00	19.80	9.25
9.5	27	166	71.00	191.50	192.00	266.60	19.90	9.80	73.00	196.50	197.00	270.60	20.20	9.80
10	28	173	73.00	195.50	196.00	274.20	20.30	10.35	75.00	200.50	201.00	278.20	20.60	10.35
11	29	179	74.00	199.50	200.00	281.80	20.70	10.90	76.00	204.50	205.00	285.80	21.00	10.90
12	30	185	76.00	203.00	203.50	289.40	21.10	11.45	78.00	208.00	208.50	293.40	21.40	11.45
12.5	31	192	78.00	206.50	207.00	297.00	21.50	12.00	80.00	211.50	212.00	301.00	21.80	12.00
13.5	32	198	79.00	210.00	210.50	304.60	21.90	12.55	81.00	215.00	215.50	308.60	22.20	12.55
1	33	205	81.00	214.00	214.50	312.20	22.30	13.10	83.00	219.00	219.50	316.20	22.60	13.10
2	34	212	82.00	218.00	218.50	319.80	22.70	13.65	84.00	223.00	223.50	323.80	23.00	13.65
2.5	35	218	83.00	222.00	222.50	327.40	23.10	14.20	85.00	227.00	227.50	331.40	23.40	14.20
3.5	36	225	85.00	226.00	226.50	335.00	23.50	14.75	87.00	231.00	231.50	339.00	23.80	14.75
4	37	232	86.00	231.00	230.50	342.60	23.90	15.30	88.00	236.00	235.50	346.60	24.20	15.30
5	38	238	88.00	235.00	234.50	350.20	24.30	15.85	90.00	240.00	239.50	354.20	24.60	15.85
6	39	245	89.00	239.00	238.50	357.80	24.70	16.40	91.00	244.00	243.50	361.80	25.00	16.40
6.5	40	252	91.00	243.50	242.50	365.40	25.10	16.95	93.00	248.50	247.50	369.40	25.40	16.95
7.5	41	258	92.00	247.50	246.50	373.00	25.50	17.50	94.00	252.50	251.50	377.00	25.80	17.50
8	42	265	94.00	252.00	250.50	380.60	25.90	18.05	96.00	257.00	255.50	384.60	26.20	18.05
9	43	272	96.00	256.00	254.50	388.20	26.30	18.60	98.00	261.00	259.50	392.20	26.60	18.60
10	44	279	97.00	260.00	258.50	395.80	26.70	19.15	99.00	265.00	263.50	399.80	27.00	19.15
10.5	45	285	98.00	264.00	262.50	403.40	27.10	19.70	100.00	269.00	267.50	407.40	27.40	19.70
11.5	46	292	100.00	268.00	266.50	411.00	27.50	20.25	102.00	273.00	271.50	415.00	27.80	20.25

Piedro® Footwear Size and Width Conversion Chart

AFO

			Narrow			Regular			Wide		
Shoe Size		Foot Size	Joint Measurements		Instep Measurements	Joint Measurements		Instep Measurements	Joint Measurements		Instep Measurements
UK	EU	Size (mm)	Width (mm) (Fig.B)	Circ (cm) (Fig.C)	Circ (cm) (Fig.D)	Width (mm) (Fig.B)	Circ (cm) (Fig.C)	Circ (cm) (Fig.D)	Width (mm) (Fig.B)	Circ (cm) (Fig.C)	Circ (cm) (Fig.D)
4	20	119 - 124	55	15.50	21.8	62	17.00	23.0	69	18.50	24.5
5	21	125 - 131	56	16.00	22.5	63	17.50	24.0	71	19.00	25.2
5.5	22	132 - 137	57	16.50	23.3	65	18.00	24.5	72	19.50	26.0
6	23	138 - 143	59	16.85	24.0	66	18.40	25.0	73	19.85	26.8
7	24	144 - 150	60	17.20	24.7	68	18.80	25.6	74	20.20	27.6
8	25	151 - 156	62	17.60	25.4	69	19.15	26.0	76	20.60	28.3
8.5	26	157 - 162	63	18.00	26.2	70	19.50	27.5	78	21.00	29.0
9.5	27	163 - 169	64	18.40	27.1	72	19.90	28.2	80	21.30	29.8
10	28	170 - 176	66	18.70	28.0	74	20.30	29.0	81	21.80	30.6
11	29	177 - 182	68	19.10	28.5	76	20.65	29.6	82	22.15	31.2
12	30	182 - 188	69	19.50	29.0	78	21.00	30.2	83	22.50	31.8
12.5	31	189 - 195	70	19.85	29.8	79	21.50	31.0	84	22.85	32.5
13.5	32	196 - 201	71	20.20	30.6	80	22.00	31.8	85	23.20	33.2
1	33	202 - 208	72	20.40	31.2	82	22.25	32.5	86	23.60	33.9
2	34	209 - 215	74	20.90	32.0	84	22.50	33.3	87	24.00	34.6
2.5	35	216 - 221	75	21.20	32.6	85	22.85	33.9	88	24.50	35.3
3.5	36	222 - 228	76	21.70	33.3	86	23.20	34.5	89	25.00	36.0
4	37	229 - 235	78	22.20	34.1	87	23.60	35.2	90	25.30	36.7
5	38	236 - 241	79	22.60	35.0	88	24.00	36.0	91	25.60	37.5
6	39	242 - 248	80	23.10	35.8	89	24.50	36.9	92	25.95	38.2
6.5	40	249 - 255	82	23.50	36.7	90	25.00	37.8	93	26.50	39.0
7.5	41	256 - 261	84	23.85	37.4	91	25.50	38.4	94	26.85	39.9
8	42	262 - 268	85	24.20	38.0	92	26.00	39.0	95	27.20	40.8
9	43	269 - 275	86	24.70	38.8	93	26.50	39.7	96	27.50	41.4
10	44	276 - 282	88	25.10	39.6	94	26.85	40.5	97	27.90	42.0
10.5	45	283 - 288	89	25.50	40.2	96	27.20	41.4	98	28.30	42.7
11.5	46	289 - 295	90	25.85	41.0	97	27.60	42.2	99	28.80	43.4




Fig A.




Fig B.




Fig C.




Fig D.

Fig A. Joint circumference

Fig B. Joint Width measurement

Fig C. Size measurement

Fig D. Short Heel circumference

Footwear Adaptions and Repairs

Piedro® footwear adaptations are carried out by our first-class Footwear Technicians using our original Piedro® high quality materials and components for the best possible match and finish.



B023 Welt Protector / Bumper



B024 External Toe Cap



B024 Sneath Heel

Code	Adaption
B001	Outside raise up to 25mm WITH a Rand
B002	Each subsequent 25mm or part thereof
B003	Outside raise up to 25mm WITHOUT a Rand
B004	Each subsequent 25mm or part thereof
B005	Heel raise up to 25mm
B006	Heel lowering
B007	Heel Seat or Valgus / Varus Wedge
B008	Stirrup, Heel, Toe Raising or Caliper Retaining Strap
B009	Additional padding including sheepskin lined tongues
B010	Lowering the upper beneath the ankle
B011	Addition of a cuff to increase the upper height
B012	Extending the opening down to the toe
B013	Conversion to Velcro® fastening per strap
B014	Addition of lambs wool lining or similar quarters
B015	Addition of high or elongated stiffener to the medial or lateral aspect
B016	Addition of a Metatarsal Bar
B017	Rocker Sole / Sole Raise
B018	Fitted Sole Wedge
B019	Floated Sole
B020	Sole Strip or Plate in composite material
B021	Sole Strip or Plate in non-composite material or Weighted Heels
B022	Rubber or steel tip to the sole or heel
B023	Welt Protector / Bumper
B024	External Protective Toe Cap / Sneath Heel
B025	Full or Half Elongated, Oblique, Floated Heel, Heel Wedge or Torque Heel

All Footwear Adaptations apply to single items of footwear unless otherwise indicated.

Footwear Adaptions and Repairs



B008 Navicular Strap



B003 Through Shoe Raise



B003 Tapered Shoe Raise

Code	Adaption
B030	Heel Socket WITHOUT Front or Back Stops
B031	Heel Socket WITH Front or Back Stops
B033	Padded T or Y Strap
B034	Double Decker T or Y Strap
B040A	A Full or 3/4 Size non-moulded Inlay in non-shock absorbent material Bontex, Plastazote etc including sheepskin Insoles
B041	Metatarsal or Valgus Support Heel Pad fitted to B040A
B042	Metatarsal or Valgus support or Heel Pad fixed to the shoe
B043	Heel Elevator, cork HDP or EVA with suitable covering
B054A	Functional foot orthosis moulded to cast-composite material
B054B	Functional foot orthosis moulded to cast-high temperature Thermoplastic
B054C	Functional foot orthosis or total contact Inlay moulded to cast
B056	Functional foot orthotic Posting
B057	Stretching of footwear

Repairs

Code	Repair
C001	Renew inside raise up to 25mm
C002	Each subsequent 25mm or part thereof
C011	Re-stitch upper per 25mm or part thereof
C013	Replace Eyelets, Ski-Hooks or D-Rings - per pair
C021	Soling of heeling in synthetic material or rubber - per pair
C025	Through soiling, synthetic, rubber or leather
C027	Renew stick on soles - per pair

All Footwear Adaptations apply to single items of footwear unless otherwise indicated.



Place your order

You can send us your order by phone, email, fax, post or online!

Please use the Order Form on page 85

Please send to:-



0845 617 1976



0845 617 1977



www.piedro-uk.co.uk



customerservices@piedro-uk.co.uk



Piedro Ltd

Unit 6 Foxwood Ind Park | Foxwood Close
Chesterfield | Derbyshire | S41 9RB
United Kingdom

Piedro®

Firststeps

Made to Order Collection

Sizes 15 - 21

Features

- Follows the natural shape of the foot
- Provides all the heel and toe space needed to allow the foot to develop naturally
- Provides stability
- Helps improve balance
- Lightweight
- Flexible, non-slip sole
- Opens up to the toe
- Made from high quality leather and components for breathability, comfort and longevity
- Minimal seams to prevent rubbing
- Available in sizes 15 - 21
- Available in widths 3½ and 5½

Firststeps Sizes 15 - 21

Sparkle

2299.9900

White Leather

Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes



Shine

2299.5600

Dark Blue Leather

Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes



Twinkle

2299.0100

Fuchsia Leather

Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes



Dazzle

2299.6500

Red Leather

Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes



Firststeps Sizes 15 - 21

Persian

2309.4537

Purple Leather

Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes



Tabby

2309.0137

Fuchsia Leather

Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes



Manx

2309.9837

Patent Black Leather

Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes



Firststeps Sizes 15 - 21

Westie

2301.5736

Cobalt Combi Leather



Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes

Sheltie

2301.6536

Red Combi Leather



Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes

Russell

2301.5636

Dark Blue Combi Leather



Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes

Carter

2301.1536

Dark Brown Combi Leather



Size	15 - 21
Widths	3½, 5½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	3½, 5½
Stability	3½, 5½
Lace	Yes
Velcro	Yes

Piedro®

Stitchdown

Made to Order Collection

Sizes 18 - 42

Features

- Extra heel and sole rigidity for increased foot support
- Flexible forefoot for natural gait development
- Supports the Malleoli and helps to control the ankle joint
- Helps to promote good posture
- Provides all the heel and toe space needed to allow the foot to develop naturally
- Contoured insole, positions the foot correctly inside the shoe
- Made from high quality leather and components for breathability, comfort and longevity
- Minimal seams to prevent rubbing
- Available in sizes 18 - 42
- Available in widths 1½ - 8½

Sole options for Stitchdown Shoes

Rubber Sole Cleo



Black

Amber

Blue

Pink

EVA Port Flex Lightweight Sole



Black

Amber

White

Stitchdown Sizes 18 - 26

Pepper

2307.5600

Dark Blue Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Astra

2307.9805

Patent Black Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Chilli

2307.6536

Red Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Anise

2307.0126

Fuchsia Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 18 - 26

Jack

2303.5680

Dark Blue Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Step

2303.9885

Black Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stride

2303.8098

Grey Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Leap

2303.1556

Dark Brown Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 18 - 26

Hop

2305.5636

Dark Blue Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Bounce

2305.1536

Dark Brown Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Pace

2305.9836

Black Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Walk

2305.8026

Grey Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 18 - 26

All models made with hard wearing heel and protective toe cap

Bumper

2340.0126

Fuchsia Nubuck



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Shimmer

2340.9805

Black Patent



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Evan

2342.9885

Black Oiled Nubuck



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Rock

2342.1536

Dark Brown Oiled Nubuck



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Stitchdown Sizes 18 - 26

Flash

2510.5637

Dark Blue Leather



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Cosmic

2510.4537

Purple Leather



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Comet

2510.6537

Red Leather



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Star

2510.0136

Fuchsia Leather



Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Stitchdown Sizes 18 - 26

Blossom

2300.5637
Black Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Flower

2300.4537
Purple Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Petal

2300.6537
Red Leather

Size	18 - 26
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes





Place your order

You can send us your order by phone, email, fax, post or online!

Please use the Order Form on page 85

Please send to:-



0845 617 1976



0845 617 1977



www.piedro-uk.co.uk



customerservices@piedro-uk.co.uk



Piedro Ltd

Unit 6 Foxwood Ind Park | Foxwood Close
Chesterfield | Derbyshire | S41 9RB
United Kingdom



Stitchdown Sizes 24 - 42

Spin

2482.9826

Black Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Swirl

2482.6536

Red Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Dart

2482.4526

Purple Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Dash

2482.0126

Fuchsia Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Stitchdown Sizes 24 - 42



Leopard

2494.5637

Dark Blue Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Cougar

2494.6537

Red Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Lynx

2494.4537

Purple Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 24 - 42

Holly

2492.5637

Dark Blue Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Crocus

2492.4537

Purple Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Daisy

2492.6537

Red Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Lily

2492.0136

Fuchsia Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Stitchdown Sizes 24 - 42



Frankie

2500.5636

Dark Blue Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Jamie

2500.9836

Black Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stevie

2500.1526

Brown Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 24 - 42

Fox

2483.9802
Black Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Wolf

2483.9800
Black Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Bear

2483.2002
Camel Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Bison

2483.1402
Dark Brown Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Deer

2483.0102
Fuchsia Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 24 - 42



Bobbie

2502.9836
Black Combi Leather

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Charlie

2502.1526
Dark Brown Combi Leather

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Chris

2502.5636
Dark Blue Combi Leather

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 24 - 42

Trek

2484.9802

Black Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Rambler

2484.9804

Black Oiled Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Trail

2484.1502

Dark Brown Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

March

2484.1501

Dark Brown Oiled Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Troop

2484.2002

Camel Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Hiker

2484.5602

Dark Blue Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Stitchdown Sizes 24 - 42

Button

2504.8025

Grey Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Verstappen

2504.5665

Dark Blue Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Hamilton

2504.1509

Dark Brown Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 24 - 42

Super

2496.9826

Black Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Bat

2496.5636

Dark Blue Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Wonder

2496.0126

Fuchsia Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Vixen

2496.4536

Purple Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Wasp

2496.6536

Red Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 24 - 42

Phantom

2498.9826

Black Combi Leather

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Grayson

2498.8026

Grey Combi Leather

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Thor

2498.1526

Dark Brown Combi Leather

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 24 - 42

Pangolin

2506.9836

Black Combi Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Mongoose

2506.8026

Grey Combi Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Jackal

2506.5636

Dark Blue Combi Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Meerkat

2506.1526

Dark Brown Combi Leather Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Stitchdown Sizes 24 - 42

Rhino

2508.9885

Black Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Giraffe

2508.1556

Dark Brown Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Elephant

2508.8056

Grey Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Buffalo

2508.2000

Cognac Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Hippo

2508.5601

Dark Blue Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 18 - 46

George

2481.9836

Black Leather Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Bryn

2481.5602

Dark Blue Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Abbie

2481.0102

Fuchsia Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Gail

2481.4502

Purple Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 18 - 46

Jeff

2480.5636

Dark Blue Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Gary

2480.5611

Dark Blue Leather

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Rebecca

2480.0126

Fuchsia Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Isabel

2480.0102

Fuchsia Leather

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



Stitchdown Sizes 18 - 46

Freddy

2480.1462

Dark Brown Nubuck Leather



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Reggie

2480.1426

Dark Brown Leather Nubuck



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Tudor

2480.9800

Black Leather



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes

Elizabeth

2480.9805

Black Patent



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Lace	Yes
Velcro	Yes



PIEDRO
always one step ahead

Marina

2485.9805

Black Patent Leather

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Velcro	Yes



Lucy

2476.9837

Black Leather Crocodile Detail

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight Sole Option	
Velcro	Yes



Place your order

You can send us your order by phone, email, fax, post or online!

Please use the Order Form on page 85

Please send to:-



0845 617 1976



0845 617 1977



www.piedro-uk.co.uk



customerservices@piedro-uk.co.uk



Piedro Ltd

Unit 6 Foxwood Ind Park | Foxwood Close
Chesterfield | Derbyshire | S41 9RB
United Kingdom

Piedro® Trainers

Made to Order Collection

Sizes 24 - 42

Features

- Extra heel and sole rigidity for increased foot support
- Flexible forefoot for natural gait development
- Supports the Malleoli and helps to control the ankle joint
- Helps to promote good posture
- Provides all the heel and toe space needed to allow the foot to develop naturally
- Contoured insole, positions the foot correctly inside the shoe
- Made from high quality leather and components for breathability, comfort and longevity
- Minimal seams to prevent rubbing
- Available in sizes 24 - 42
- Available in widths 1½ - 8½

Sole options for Trainers

Rubber Outsole Runner



Black

Amber

Available in EVA wedge colors:



Black



Grey



White

Trainers Sizes 24 - 42

Graff

2034.9957

White Cobalt Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Moore

2034.9960

White Blue Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	EVA
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Lace	Yes
Velcro	Yes

Coe

2034.5736

Cobalt Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Ennis

2034.0336

Fuchsia Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Lennox

2034.9885

Black Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Trainers Sizes 24 - 42

Gunnell

2036.0126

Fuchsia Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Grainger

2036.4526

Purple Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Cram

2036.5736

Cobalt Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Foster

2036.6557

Red Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Wiggins

2036.9885

Black Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Drew

2036.9966

White Red Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Trainers Sizes 24 - 42

Nadia

2038.4526

Purple Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Novak

2038.9836

Black Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Melissa

2038.0126

Fuchsia Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Lionel

2038.5636

Dark Blue Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Rory

2038.9936

White Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Trainers Sizes 24 - 42

Emma

2027.4526

Purple Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Vettel

2027.9885

Black Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Venus

2027.0126

Fuchsia Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Neymar

2027.5736

Cobalt Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Alonso

2027.9957

White Combi Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes





PIEDRO
always one step ahead



Piedro[®]

High Tops & Mid Tops

Made to Order Collection

Sizes 24 - 42

Features

- Extra heel and sole rigidity for increased foot support
- Flexible forefoot for natural gait development
- Supports the Malleoli and helps to control the ankle joint
- Helps to promote good posture
- Provides all the heel and toe space needed to allow the foot to develop naturally
- Contoured insole, positions the foot correctly inside the shoe
- Made from high quality leather and components for breathability, comfort and longevity
- Minimal seams to prevent rubbing
- Available in sizes 24 - 42
- Available in widths 2½ - 8½

Sole options for Trainers

Combined Sole and PU Bumper
PU Bumper outsole



Black

Grey

Rubber Sole
Fish



Black

Amber

Extra durable Rubbersole
Tyres



Black

Amber

EVA Sole - **Astro Start**
Lightweight



Black

Amber

High Tops & Mid Tops Sizes 24 - 42

Limit

2124.5700

Cobalt Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Lennon

2124.9800

Black Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Martina

2124.0100

Fuchsia Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Murray

2124.1500

Dark Brown Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

High Tops & Mid Tops Sizes 24 - 42

Starr

2105.9800

Black Leather

Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Federer

2105.5612

Dark Blue Leather

Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Mini

2105.0112

Fuchsia Leather

Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Bentley

2105.1512

Dark Brown Leather

Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



High Tops & Mid Tops Sizes 24 - 42

Bolt

2123.5736

Cobalt Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Emily

2123.0126

Fuchsia Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Quartz

2123.9807

Black Yellow Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Sue

2123.9826

Black Fuchsia Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Linford

2123.1900

Camel Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Ringo

2123.9800

Black Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

High Tops & Mid Tops Sizes 24 - 42

McCartney

2114.9800

Black Leather
Velcro Model



Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Axis

2114.5700

Cobalt Leather
Velcro Model



Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Serena

2114.0100

Fuchsia Leather
Velcro Model



Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Samprass

2114.1500

Dark Brown Leather
Velcro Model



Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	Yes
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

High Tops & Mid Tops Sizes 24 - 42

Redgrave

2133.5736

Cobalt Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Anne

2133.0126

Cobalt Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Johnson

2133.9807

Black Yellow Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Jewel

2133.9826

Black Fuchsia Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Lineker

2133.1900

Camel Combi Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

Courtney

2133.9800

Black Leather



Size	24 - 42
Width	2½ - 8½
Sole Unit	Rubber
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes

High Tops & Mid Tops Sizes 24 - 42

Chevvy

2115.5612

Dark Blue Leather
Velcro Model

Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Harrison

2115.9800

Black Leather

Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Dodge

2115.0112

Fuchsia Leather
Velcro Model

Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes



Jeep

2115.1512

Dark Brown Leather
Velcro Model

Size	24 - 42
Width	2½ - 8½
Sole Unit	EVA
Sole Change	No
Rehabilitation	2½, 4½, 6½, 8½
Stability	2½, 4½, 6½, 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lace	Yes
Velcro	Yes





Place your order

You can send us your order by phone, email, fax, post or online!

Please use the Order Form on page 85

Please send to:-



0845 617 1976



0845 617 1977



www.piedro-uk.co.uk



customerservices@piedro-uk.co.uk



Piedro Ltd

Unit 6 Foxwood Ind Park | Foxwood Close
Chesterfield | Derbyshire | S41 9RB
United Kingdom

Piedro® Sandals

Made to Order Collection

Sizes 24 - 42

Features

- Extra heel and sole rigidity for increased foot support
- Flexible forefoot for natural gait development
- Supports the Malleoli and helps to control the ankle joint
- Helps to promote good posture
- Provides all the heel and toe space needed to allow the foot to develop naturally
- Contoured insole, positions the foot correctly inside the shoe
- Made from high quality leather and components for breathability, comfort and longevity
- Minimal seams to prevent rubbing
- Available in sizes 24 - 42
- Available in widths 1½ - 8½

Sole options for Sandals

Rubber Sole Cleo



EVA-Outsole - Port Flex Lightweight Sole



Sandals Sizes 24 - 42



Ibiza

2606.6536

Red Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Majorca

2606.0126

Fuchsia Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Menorca

2606.4536

Purple Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Sandals Sizes 24 - 42

Jamaica

2608.9836

Black Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Aruba

2608.1536

Dark Brown Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Sumatra

2608.8036

Grey Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Trinidad

2608.5636

Dark Blue Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Sandals Sizes 24 - 42

Bonaire

2610.9836

Black Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Kitts

2610.1536

Cognac Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Barbados

2610.8036

Grey Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Curacao

2610.5636

Dark Blue Combi Leather Nubuck

Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Sandals Sizes 24 - 42

Culebra

2514.5602

Dark Blue Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Antigua

2514.0136

Fuchsia Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Nevis

2514.9802

Black Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Tortola

2514.1502

Dark Brown Nubuck



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Belize

2514.4537

Purple Leather



Size	24 - 42
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Sandals Sizes 18 - 46



Sandy

2601.1512

Dark Brown Oiled Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Centurion

2601.5680

Dark Blue Leather Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Pebble

2601.1202

Beige Leather Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Sandals Sizes 18 - 46



Wave

2602.5680

Dark Blue Leather Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Java

2602.0136

Fuchsia Leather Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Nemo

2602.6537

Red Leather Nubuck

Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Sandals Sizes 18 - 46

Shingle

2604.5673

Dark Blue Combi Nubuck Leather



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Coast

2604.1509

Dark Brown Leather



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Summer

2604.0136

Fuchsia Nubuck Leather



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Shelley

2604.6537

Red Nubuck Leather



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes

Palm

2604.1261

Beige Mint Combi Nubuck Leather



Size	18 - 46
Width	1½ - 8½
Sole Unit	Rubber
Sole Change	Yes
Rehabilitation	1½ - 8½
Stability	1½ - 8½
AFO	N-R-W
Reverse Lasted	4½, 6½
Lightweight	Yes
Buckle	Yes
Velcro	Yes



Piedro® Only AFO Made to Order Collection

Features

- Made on specialist footwear lasts to ensure the correct position of the lower limb is at 90°
- Additional space provided at the front and around heel
- Can accommodate stock or custom Ankle Foot Orthosis (AFO)
- Adaptable - See pages 17 - 18 for footwear adaptations
- Can accommodate an inner-shoe such as Rainier Boot or Arizona splint for example

AFO Backflips Sizes 24 - 36

Christiano

2760.9957

White Cobalt Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Colour Options



2760.0336



2760.5736



2760.0885



2760.9960



Velcro Model

AFO Backflips Sizes 24 - 36

Beckham

2761.5736

Cobalt Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Colour Options



2761.0126



2761.4526



2761.9885



2761.9957



Velcro Model

AFO Backflips Sizes 24 - 36

Cruyff

2762.9957

Dark Blue Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Colour Options



2762.0126



2762.4526



2762.9936



2762.9836



Velcro Model

AFO Shoes Sizes 24 - 36

Anya

2138.0126

Fuchsia Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Heather

2138.4526

Purple Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Bairstow

2138.9936

White Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Geraint

2138.5636

Blue Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Moeen

2138.9836

Black Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



AFO Shoes Sizes 24 - 36

Williams

2134.0336

Fuchsia Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Cash

2134.5736

Cobalt Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Lendl

2134.9960

White Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Edberg

2134.9957

White Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Hewitt

2134.9885

Black Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



AFO Shoes Sizes 24 - 36

Henin

2127.0126

Fuchsia Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Agassi

2127.4526

Purple Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Becker

2127.0126

Fuchsia Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Hingis

2127.5736

Cobalt Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes



Sanchez

2127.9957

White Combi Leather

AFO	Widths: N - R - W
Size	24 - 36
Sole Unit	Rubber
Sole Change	No
Velcro	Yes
Lace	Yes





PIEDRO
always one step ahead

Place your order

You can send us your order by phone, email, fax, post or online!

Please use the Order Form on page 85

Please send to:-



0845 617 1976



0845 617 1977



www.piedro-uk.co.uk



customerservices@piedro-uk.co.uk



Piedro Ltd

Unit 6 Foxwood Ind Park | Foxwood Close
Chesterfield | Derbyshire | S41 9RB
United Kingdom



Insoles Accessories Footwear Care




Insoles

Piedro® produces a range of off-the-shelf and custom insoles .

Our RESPRO® and BioFit® Inlays are made on the same lasts as our Piedro® footwear, ensuring a perfect fit every time.

Normal Arch



Insole/Inlay

- BioFit Pronation
- Heel Cups
- BioFit Neutral Inlay
- RESPRO ¾ Size
- RESPRO Full Size


Use

Prevention
Pronation Control

Footwear Type

Rehabilitation

Flat Arch




- Heel Cups
- BioFit Neutral Inlay
- RESPRO ¾ Size
- RESPRO Full Size

Mild to moderate
over Pronation
Flat foot

Rehabilitation

High Arch



- BioFit Dynamic Inlay
- RESPRO Arch Filled ¾ Size Inlay
- RESPRO Arch Filled Full Size Inlay

Moderate to severe
over Pronation
Supination

Stability

Heel Cups



Mild to moderate Over Pronation
Supplied in pairs



Features

- Available in sizes 18 – 44 and supplied in pairs
- 3cm heel cup depth
- Compatible with Pedro® and most shop bought footwear
- Provides support and stabilization for the developing foot
- Contoured to provide support to the longitudinal and peroneal arches
- Helps to support and stabilize the subtalar and midtarsal joints
- Neutral rear foot post incorporated
- Easy addition of posts, wedges or raises
- Navicular padding supplied

Item	Size	Width
18 - 26	3, 5, 7	HC01
27 - 44	3, 5, 7	HC02

BioFit® Functional Inlays

BioFit® - Biomechanical Functional Inlay Triple System

Our BioFit® Inlay range is fully compatible and recommended for use with all Pedro® footwear.

Item	Size	Width	Codes
1. Dynamic Inlay	18 - 37	1½, 2½, 3½, 4½, 5½, 6½, 7½, 8½	98001
2. Neutral Inlay	18 - 27	1½, 2½, 3½, 4½, 5½, 6½, 7½, 8½	98000
3. Pronation Control	18 - 37	1½, 2½, 3½, 4½, 5½, 6½, 7½, 8½	98002
4. Inlay Posting Set	3° & 5°	-	98003

BioFit® Dynamic Inlay

Overpronation, Cerebral Palsy, Neurological Impairment, support for longitudinal & peroneal arches

1.



Features

- Split sizes and widths available with no surcharge
- Provides support and stabilization for the developing foot
- Contoured to provide support to the longitudinal and peroneal arches without undue pressure
- Deep heel cup with medial lateral flanges
- Neutral rearfoot post incorporated into moulding
- Stabilizes and supports the tarso-metatarsal junction and the metatarsal dome
- Stabilizes the metatarsal phalangeal joints, helps to maintain the integrity of the transverse curvature
- Stabilizes the hallux, promoting better reactions
- Easy addition of posts, wedges or raises
- Fully customizable to your patients' requirements
- 3° or 5° wedging strip available – please request when ordering
- Single construction, heat mouldable EVA shell
- Suedette cover

BioFit® Functional Inlays

2.



BioFit® Neutral Inlay

Support and correction for the overpronating foot

Features:

- Split sizes and widths available with no surcharge
- Contoured to provide support to the longitudinal and peroneal arches without undue pressure
- Deep heel cup. When used in combination with Pedro® footwear helps to support and stabilize the midtarsal joints
- Easy addition of posts, wedges or raises
- Customizable to your individual patient's requirements
- 3° or 5° wedging strip available – please request when ordering
- Single construction, heat mouldable EVA shell
- Suedette cover

3.



BioFit® Pronation Control Inlay

Support and stabilization of the developing foot

Features:

- Rigid structure
- Split sizes and widths available with no surcharge
- Contoured to provide support to the longitudinal and peroneal arches
- Deep heel cup. When used in combination with Pedro® footwear helps to support and stabilize the mid- tarsal joints
- Easy addition of posts, wedges or raises
- Customizable to your individual patients requirements
- 3° or 5° wedging strip available - please request when ordering
- Heat mouldable shell Suedette cover

4.



BioFit® Neutral Inlay

Support and correction of the overpronating foot

Features:

- Rigid structure
- Split sizes and widths available with no surcharge
- Contoured to provide support to the longitudinal and peroneal arches
- Deep heel cup. When used in combination with Pedro® footwear helps to support and stabilize the mid- tarsal joints
- Easy addition of posts, wedges or raises
- Customizable to your individual patients requirements
- 3° or 5° wedging strip available - please request when ordering
- Heat mouldable shell Suedette cover

Respro® Standard and Arch Filled Inlays

Our most popular inlay range, available in full and $\frac{3}{4}$ sizes with standard and arch filled options, which provide support for the over-pronating foot.

Item	Size	Width	Codes
1. Standard Full Size	18 - 50	1½, 2½, 3½, 4½, 5½, 6½, 7½, 8½	98004
2. Standard $\frac{3}{4}$ Size	18 - 45	1½, 2½, 3½, 4½, 5½, 6½, 7½, 8½	98005
3. Arch Filled Size	18 - 50	1½, 2½, 3½, 4½, 5½, 6½, 7½, 8½	98006
4. Arch Filled $\frac{3}{4}$ Size	18 - 45	1½, 2½, 3½, 4½, 5½, 6½, 7½, 8½	98007

RESPRO® - Standard full and $\frac{3}{4}$ size inlays



RESPRO® - Arch filled full and $\frac{3}{4}$ size inlays



RESPRO® Standard Inlays

Mild to Moderate Overpronation

Features

- Full size and $\frac{3}{4}$ size available
- Standard and Arch Filled options
- Compatible with Pedro® and most shop bought footwear
- Contoured to provide support to the longitudinal and peroneal arches
- Neutral rearfoot post incorporated into moulding
- Easy addition of posts, wedges or raises
- Stabilizes the metatarsal phalangeal joints, helps to maintain the integrity of the transverse curvature
- Deep heel cup, when used with Pedro® footwear helps to stabilize the subtalar and midtarsal joints
- Easy addition of posts, wedges or raises
- Fully customizable to your patients' requirements
- 3° or 5° wedging strips available
- please request when ordering
- Heat mouldable

RESPRO® Arch Filled Inlays

Moderate to High Levels of Overpronation

Arch Filled RESPRO® Inlays provide a higher level of arch support than our standard model, as the arch area is filled.

Arch Filled shoes are recommended where higher levels of foot control and support are required.

Features:

- Holds to relieve foot, ankle, leg, hip and lower back pain related to overpronation
- RESPRO® inlays fit perfectly into Pedro® footwear
- Holds foot in a corrected anatomical position
- Helps to distribute weight evenly over the feet
- Available in full and $\frac{3}{4}$ Size models

Accessories

Take good care of your Pedro® footwear. Here is a selection of our recommended footwear care products.



Code: TC001
Collonil Clean & Care
textile cleaner



Code: 99001
Pedro® Foot Measure

Code: WS001
Collonil Waterstop
shoe cream



Code: CR001
Collonil Cleaning Rubber



Code: 99070
Replacement Pedro®
shoe laces



Pedro® Insoles

Widths	Size	Code
1½ - 2½	18 - 46	99050
4½ - 5½	18 - 46	99055
6½ - 7½	18 - 46	99060



Shoe polish

Black	PFC01
Dark Blue	PFC02
Red	PFC03

Please state Model, Size and Width required when ordering.

Order Number Specialist

Address

Client Last Details Date




Piedro® Concept Range - Made to Order
PMF® - Piedro® Modular Footwear

Model & Sizing

Model Number Colour Number Size Width Stability Rehabilitation AFO Reverse Lasted

Feet & Lower Leg Measurements



	Left	Right
1. Foot Size		
2. Joint Width		
3. Joint Circumference		
4. Instep Circumference		
5. Short Heel Circumference		
6. Toe Box Height		


Custom Leg and Ankle Circumference (Circumference Adaption)

Left	Height	Right
	10cm	
	15cm	
	20cm	
	25cm	

Fitting Stage Required ☐

Additions

Last



	Left	Right
1. Lateral Joint Width		
2. Medial Joint Width		
3. Lateral Heel Width		
4. Medial Heel Width		
5. Hammer Toe		
6. Bunionette		

	Left	Right
7. Hallux		
8. Depth to Forepart		
9. Depth to Toe Heel		
10. Depth Over Instep		
11. Straighten Heel Clip		
12. Heel Only		

Upper Adaptions

Lining

Leather Synthetic Fur Real Fur Diabetic Sympatex

Closure Laces

Eyelets Hooks D-Rings Blind Eyelets Buckle & Strap

Closure Velcro & Straps

Return Velcro Retaining Strap Lap-Over Velcro Navicular/U Strap Single Hand Velcro Caliper Retaining Strap Velcro Separate

Toe Puffs

Soft Hard Rim

Heel Stiffener

Stability Rehabilitation Heel Exostosis None


Upper Straps

T Strap Y Strap Double Decker Strap

Stiffener Hardness

Soft Standard Hard Double Extra Padding

Stiffener Length



Allergy

Lining Upper Metal Free

Adult Footwear Only

Stretch Leather

Sole & Heel Adaptions

Custom Raise

	Left	Right
Toes		
Joint		
Heel		

Sole Float

Left Medial	Right Medial
Left Lateral	Right Lateral

Heel Float

Left Medial	Right Medial
Left Lateral	Right Lateral

Sole Wedge

Left Medial	Right Medial
Left Lateral	Right Lateral

Heel Wedge

Left Medial	Right Medial
Left Lateral	Right Lateral

Rocker Sole Types

Normal Rocker




Advancing Rocker



Polyphase Rocker



Delaying Rocker



2-Phase Rocker



Removable Carbon Insole Full Carbon Fibre Sole Plate Carbon Sole Strip Separate Outsoles

Sneath Heel Buffer Heel Sach Heel Thomas Heel

Heel Socket

	Medial	Lateral	Double		Medial	Lateral	Double
Tube 5/16"				Tube 9/16" x 3/16"			
Tube 1/4"				Tube 1/2" x 1/8"			

Other

Welt Protector Protective Toe Cap Elongation Medial Lateral Extra Pair of Laces Fast Delivery (5 Working Days)

Additional Comments

Order Information

Efficient delivery

We aim to dispatch all stock item orders on the day of receipt. If you place your order before 3 p.m. on Monday to Friday, we will always try to ship your order on the same day, or otherwise the next working day. Made to Order or Out of Stock items will take longer to dispatch. Due to our wide variety of footwear styles, in numerous sizes, and colour and width options, we occasionally encounter stock shortages on popular items. In this case, we will endeavour to offer an alternative item, or otherwise dispatch your order within 2 to 3 weeks.

All orders are shipped according to the customer's request but this cannot be guaranteed. From time to time, orders may be delayed due to external factors, such as the carrier or the weather.

Payment terms

All payments are in GB Pounds (£). Payment is due on ordering the items, prior to dispatch. For customers with a Piedro® Credit Account, payment is due within 30 days of the date of our invoice. Please see our price list for our full payment terms and conditions.

Making a complaint

If you wish to make a complaint, please get in touch with our Customer Services Department first. They will be glad to assist you with any queries or problems.

Returns Policy

If for any reason you are dissatisfied with your items, simply complete the returns form (*on the reverse of your delivery note*) and send your unused/unworn items to us in their original packaging within 30 days of the order date for a refund or exchange. Unfortunately any footwear or insoles which have been made to your specific requirements is considered as non-refundable, unless it is considered faulty. Please see our price list for our full Returns Policy.

Made to Order, Split Size and Width Footwear

All made to order items, including split size and width footwear, is non-refundable, unless faulty. Dispatch of made to order items is usually within 2 to 3 weeks of the order date.

Extra Catalogues and Price Lists

If you or any of your friends and colleagues would like a copy of our catalogue or price list, please call our Customer Services Department on **0845 617 1976** or email **customerservices@piedro-uk.co.uk**

Additional information

From time to time, styles, colours and prices may vary slightly. Slight colour variations may also occur during the printing process.

Delivery Options

For detailed information on our delivery services, rates and conditions, please contact our Customer Services Department.



0845 617 1976



0845 617 1977



www.piedro-uk.co.uk



customerservices@piedro-uk.co.uk



Piedro Ltd

Unit 6 Foxwood Ind Park | Foxwood Close
Chesterfield | Derbyshire | S41 9RB
United Kingdom

Piedro[®]

is one of the world's leading
manufacturers of high-quality
Therapy & Orthopaedic
Footwear.

Order online, phone, post or fax



p: 0845 617 1976
f: 0845 617 1977
www.piedro-uk.co.uk
customerservices@piedro-uk.co.uk

Piedro Ltd

Unit 6 Foxwood Ind Park | Foxwood Close
Chesterfield | Derbyshire | S41 9RB
United Kingdom



PIEDRO®
alwaysonestepahead